



Inno⁺Fine

Biopsy Needle Guides for MINDRAY

- Integral design without any accessories
- High grade medical stainless steel
- Artificial polishing
- Supporting accurate and angle adjustable biopsy
- Matching a wide range of needle gauges
- Simple pressing, quick release

List of Biopsy Needle Guides for MINDRAY

REF#	Transducer	Tranducer Type	Gauge Size/Guide Channel Length	REF#	Transducer	Tranducer Type	Gauge Size/Guide Channel Length
JSM-015	C5-2/C5-2s/C5-2E	Convex	11-23G,3Fr-9Fr/3.5cm	JSM-127	L14-5W/L14-5WU(ZONARE)	Linear	11-23G,3Fr-9Fr/3.5cm
JSM-006	3C5A/3C5/3C5s/3C5P	Convex	14-22G,3Fr-6Fr/3.5cm	JSM-255	L8-3(ZONARE)	Linear	14-22G,3Fr-6Fr/3.5cm
JSM-001	35C50EA/35C50EB/35C50P	Convex	11-23G,3Fr-9Fr/3.5cm	JSM-264	L9-3U	Linear	11-23G,3Fr-9Fr/3.5cm
JSM-045	SC5-1E	Convex	11-23G,3Fr-9Fr/3.5cm	JSM-074	3C1/3C1s	Microconvex	14-22G,3Fr-6Fr/3.5cm
JSM-022	C5-1s/C5-1U/SC6-1U	Convex	11-23G,3Fr-9Fr/3.5cm	JSM-075	65C15EA/6C2/6C2s/6C2P	Microconvex	16/18G,3.2cm
JSM-219	C7-3E	Convex	11-23G,3Fr-9Fr/3.5cm	JSM-118	C11-3s/C11-3E/C11-3U	Microconvex	11-23G,3Fr-9Fr/3.5cm
JSM-031	SC5-1U	Convex	11-23G,3Fr-9Fr/3.5cm	JSM-024	C6-2GE/C6-2Gs/C6-2GU	Microconvex	14-22G,3Fr-9Fr/3.5cm
JSM-116	C6-2(ZONARE)	Convex	16/18G,3.2cm	JSM-123	C4-1(ZONARE)	Microconvex	16/18G,3cm
JSM-287	L20-5/L20-5s/L20-5U	Linear	11-23G,3Fr-9Fr/3.5cm	JSM-229	C9-3(ZONARE)	Microconvex	16/18G,3.2cm
JSM-036	LM15-6U	Linear	11-23G,3Fr-9Fr/3.5cm	JSM-265	SC8-2U	Microconvex	11-23G,3Fr-9Fr/3.5cm
JSM-146	L12-3RC	Linear	20-23G/1cm	JSM-011	2P2(s)/2P2P/P4-2(s)/P4-2E/P4-2NE/SP5-1s/SP5-1E/SP5-1U	Phased Array	14-22G,3Fr-6Fr/3.5cm
JSM-284	L10-3E/L10-3s	Linear	14-22G,3Fr-6Fr/3.5cm	JSM-221	P7-3/P7-3s/P7-3E	Phased Array	11-23G,3Fr-9Fr/3.5cm
JSM-214	L14-3U	Linear	11-23G,3Fr-9Fr/3.5cm	JSM-144	D8-4U	4D	16/18G,3.2cm
JSM-007	7L4(s)/7L4A/7L4P/10L4(s)/7L5(s)/7L5P/L7-3(s)/L7-3E/L11-4(s)/L12-4(s)/L12-3E/L14-6Ns/L14-6NE/L14-6Ws/L14-6WE/L14-6WU/75L53EA(HA)	Linear	14-22G,3Fr-6Fr/3.5cm	JEM-028	65EB10EA/65EC10EA/65EC10EB/65EC10EC/65EC10ED/6CV1(s)/6CV1P/CB10-4(s)/CB10-4P/CB10-4E/CB10-4A/V10-4(s)/V10-4B(s)/V10-4BP/V11-3Ws/V11-3(E)/V11-3B(E)/V11-3W(E)/E9-3/V11-3HU/V11-3HE	Transvaginal	16-18G/14.5cm
JSM-088	75L38EA/75L38EB/75L38P	Linear	11-23G,3Fr-9Fr/3.5cm	JEM-042	E9-4(ZONARE)	Transvaginal	16-18G/14.5cm
JSM-237	LM14-6E/LM14-6s/LM16-4U	Linear	11-23G,3Fr-9Fr/3.5cm	JEM-063	DE10-3/DE10-3E/DE10-3U/DE10-3WU	Transvaginal 4D	16-18G/19cm
JSM-274	L14-6(s)/L14-6P/L14-6E/10L24EA	Linear	14-22G,3Fr-6Fr/3.5cm	JSM-109	6LB7(s)/6LB7P/6LE7(s)/6LE7P/6SEL60EA	Biplane	14/16/18/20G,3cm
JSM-126	L11-3U	Linear	11-23G,3Fr-9Fr/3.5cm				
JSM-204	L16-4Hs/L16-4HE	Linear	16/18G,3.2cm				

(Please contact InnoFine for more models.)