



InnoFine

# Biopsy Needle Guides for SONOSCAPE



- Integral design without any accessories
  - High grade medical stainless steel
  - Artificial polishing
- Supporting accurate and angle adjustable biopsy
- Matching a wide range of needle gauges
  - Simple pressing, quick release

## List of Biopsy Needle Guides for SonoScape

REF#	Transducer	Tranducer Type	Gauge Size/Guide Channel Length
JSM-183	C351	Convex	11-23G,3Fr-9Fr/3.5cm
JSM-185	C344	Convex	11-23G,3Fr-9Fr/3.5cm
JSM-190	3C-A	Convex	11-23G,3Fr-9Fr/3.5cm
JSM-206	C1-6	Convex	11-23G,3Fr-9Fr/3.5cm
JSM-208	C353	Convex	11-23G,3Fr-9Fr/3.5cm
JSM-306	C1-6A	Convex	14-22G,3Fr-6Fr/3.5cm
JSM-277	10L-I	Linear	11-23G,3Fr-9Fr/3.5cm
JSM-216	12L-A	Linear	11-23G,3Fr-9Fr/3.5cm
JSM-275	12L-B	Linear	11-23G,3Fr-9Fr/3.5cm
JSM-152	L752	Linear	11-23G,3Fr-9Fr/3.5cm
JSM-171	9L-A	Linear	11-23G,3Fr-9Fr/3.5cm
JSM-207	L742	Linear	11-23G,3Fr-9Fr/3.5cm
JSM-203	L741	Linear	11-23G,3Fr-9Fr/3.5cm
JSM-166	C322	Microconvex	11-23G,3Fr-9Fr/3.5cm

REF#	Transducer	Tranducer Type	Gauge Size/Guide Channel Length
JSM-186	C612	Microconvex	11-23G,3Fr-9Fr/3.5cm
JSM-189	C613	Microconvex	11-23G,3Fr-9Fr/3.5cm
JSM-135	4P-A	Phased Array	11-23G,3Fr-9Fr/3.5cm
JSM-205	7P-A	Phased Array	14-22G,3Fr-6Fr/3.5cm
JSM-231	S1-5	Phased Array	14-22G,3Fr-6Fr/3.5cm
JEM-017	6V1(Old version)	Transvaginal	16-18G,15cm
JEM-006	6V1(New version)	Transvaginal	16-18G,15cm
JEM-003	6V3(Old version)	Transvaginal	16-18G,15cm
JEM-031	6V3(New version)	Transvaginal	16-18G,14.5cm
JEM-009	6V4	Transvaginal	16-18G,14.5cm
JEM-047	6V7	Transvaginal	16-18G,21cm
JEM-046	VE9-5	Transvaginal	16-18G,19cm
JEM-053	BCC9-5/BCC9-4	Transvaginal	16-18G,15cm
JEM-052	BCL10-5	Transvaginal	14/16G/18/20G,3cm

(Please contact InnoFine for more models.)