

Differences in Inobitech VS OsiriX

Sr	Inobitech	OsiriX MD	Comment
1	DICOM Print with custom Layouts	Not Avalibale	With Pcare RIS also One More Print License so you get 2 licenses with us
2	Angio Structured Reporting: Vessle Analysis Module	Not Avalibale, Manual Only, One vessle at a time in curved MPR	Inobitech Scores
3	SUV Measurements: Formulae available LBM (James), LBM, BSA, BM	Only Body Mass	Inobitech Scores
4	Measurement Tools: Length, Angle, Rectangle, Oval,Open Polygon, Freehand Pencil, Point, Brush, Perpendicular Lines Angle,	Length, Angle, Rectangle, Oval,Open Polygon,Closed Polygon, Freehand Pencil, Point, Brush, Perpendicular Lines Angle, Dynamic angle,	Inobitec Don't Have Freehand Pencil & Perpendicular line ROI used for TT/TG Measurement tool, But can be done with Lasso tool
5	Inobitech Lacks feature but same ojective can be achived using camera shot tool that collects images for printing	Key Image & Single click export/Mail But uplimatealy osiriX can't print it!	Both softwares are having Pros & cons
6	Reading	Reading	Both Similar
7	Zoom/Pan/Window	Zoom/Pan/Window	Both Similar
8	Distance And Angel	Distance And Angel	Both Similar
9	Marker In Comparison	Marker In Comparison	Both Similar
10	Region Of Interest	Region Of Interest	Both Similar
11	Volume Of Interest (Freehand, Sphere)	Volume Of Interest (Freehand, Sphere)	Both Similar
12	Arrow, Pixel Lens	Arrow, Pixel Lens	Both Similar
13	Synchronized Scrolling	Synchronized Scrolling	Both Similar
14	Correlated Cursors	Correlated Cursors	Both Similar
15	MPR (Thin / Thick)	MPR (Thin / Thick)	Both Similar
16	MIP (Thin / Thick)	MIP (Thin / Thick)	Both Similar
17	MinIP	MinIP (Thin/Thick)	Both Similar OsiriX Bit scores
18	VRT (Default)	VRT (Default)	Both Similar
19	VRT User Specific Presets	Only standard Presets	Inobitech Scores
20	VRT Gallery	VRT Gallery	Both Similar
21	Clip Plane (Slice/Slab)	Cutting within clipbox	Both Similar
22	Clip Box, Punching	Clip Box, Punching	Both Similar
23	One Click Bone Removal	One Click Bone Removal	Both Similar

24	Bone Removal With Threshold	Bone Removal With Threshold Adjustment	Both Similar But methods are different Inobitech shows preview of what will be deleted
25	Single click Table Removal	Not Avalibale	Inobitech Scores
26	Parallel And Radial Ranges	Parallel And Radial Ranges	
27	Curved Ranges	Curved Ranges	
28	2D & 3D Reference Lines	2D & 3D Reference Lines	Both Similar
29	3D Reference Point	3D Reference Point	Both Similar
30	Movie (incl. export)	Movie (incl. export)	Both Similar
31	Interactive Segmentation	Interactive Segmentation	Inobitech Has More comprehensive tools
32	Region Growing / Vessel Growing	Region Growing / Vessel Growing	Both Similar
33	Non Dicom (jpeg) Saving	Non Dicom (jpeg, PDF Screenshot) Saving	OsiriX Scores
34	Manual Vessel Tracking	Manual Vessel Tracking	Both Similar
35	Case Navigator	Case Navigator	Both Similar
36	Image Substraction	Image Substraction	Both Similar
37	Color LUT in MPR	Color LUT in MPR	Both Similar
38	Image Fusion	Image Fusion	Both Similar
39	Multi Time Point Comparison	Multi Time Point Comparison	Both Similar
40	Multi Patient Comparison	Multi Patient Comparison	Both Similar
41	Multi Modality Comparison	Multi Modality Comparison	Both Similar
42	Dental	Dental	Both Similar
43	Flythrough / Bronchoscopy / Colonoscopy	Flythrough / Bronchoscopy / Colonoscopy	Both Similar
44	Measurements On Bronchoscopy Images	Measurements On Bronchoscopy Images	Both Similar
45	Labelling On Bronchoscopy Images	Labelling On Bronchoscopy Images	Both Similar
46	Spine Screw Planning Module	Not Avalibale	Inobitech Scores
47	Cobb's Angle Not Available But can be done	Cobbs Angle	OsiriX Scores But Inobitec also can do that
48	Calcium Score	Calcium Score	Both Similar
49	Anonymize cases	Anonymize cases	Both Similar
50	Series Cine Loop Play Shortcut Muple clicks	Series Cine Loop Play Shortcut Muple clicks	OsiriX Scores
51	Drag & Drop Images Directly to word Powerpoint NOT POSSIBLE, Save to Jepeg & Then Insert	Drag & Drop Images Directly to word Powerpoint directly possible	OsiriX Scores