


# ACQUIDR

DR Upgrade Solution By Retrofit

## Advanced Digital Image Acquisition System



### Simple and Efficient Solution for Your Digital Radiography

ACQUIDR is the digital imaging system composed of a Flat Panel Detector(FPD) and an imaging workstation with software. The digital FPD and fully featured imaging software with excellent digital image processing meets most of needs in diagnostic digital radiographic field.

#### Product Component

##### Wireless FPD



##### Flat Panel Detector

Imaging Area	14x17inch	17x17inch
Scintillator	CsI / GOS	
Type	Wireless	Wired

##### Imaging Software & Workstation




\*Laptop optional

## Feature

- Adaptable with most existing analog x-ray systems
- AED(Auto Exposure Detection) function available
- Simple and convenient installation
- Intuitive graphic user interface
- Over 1300 pre-loaded exam profiles with DRGEM's generator
- Configurable automatic study
- PMS / EMR / HIS / RIS work list support
- Multiple printing formats
- DICOM 3.0 compliant
- Brightness and contrast adjustment
- Rotation
- Horizontal and vertical flip
- Zoom and pan
- Left / right makers
- Optimized work flow
- Automatic shutters
- Support DAP interface
- Exposure index (EI)

### Robust and Durable Detector




Detector water & dust proof (IP51)

### High image Quality




## DR Upgrade with GXR-S

Conventional X-ray System (Film)



Digital X-ray System (DR)



©DRGEM Corporation - All rights reserved. Reproduction in whole or in part is prohibited without the written consent of the copyright-holder.

Specifications and features are subject to change without notice. All changes will be in compliance with regulations governing manufacture of medical equipment.